

## Product Description

---


The idea behind the design of the Q-Back is to allow for seamless integration of the ADI back support line to the Quantum TB3 Power Positioning System. This combination will grant the consumer an array of back support height, width, depth and profile options for their multifunctional Quantum TB3 System, resulting in a highly sophisticated seating package.



## Features

---

- The new Q-Back includes a high quality Lycra cover to complement the anti-shear features of the TB3 system, as well as offering water resistance and breathability.
- The hole pattern on the Q-Back is designed for ease of installation to the TB3 back plate as well as offering height adjustability in .5" increments equaling a total of 3.5" of height adjustment.
- Additional slots have been included on the Q-Back to allow for integration of Stealth laterals and positioning straps.
- Pre-drilled holes for mounting head-rests.
- New XX-Tall (22") height option.



*Easy flush mounting*

## Sizes

- Heights – Standard (16"), Tall (18"), X-Tall (20"), XX-Tall (22")
- Widths – 16", 18", 20"
- Contours – Active (2.5"), Deep (4"), Extra Deep (6")
- Profiles – AL Series (Scapular Cut-Out), AFT Series (Flat Top)
- **Custom options available**



Pre-drilled holes for headrest



Flush mounting



3.5" in height adjustment



High quality Lycra cover

## Coding

- **Active (2.5") Contour**  
**E2613** - POSITIONING WHEELCHAIR BACK CUSHION, POSTERIOR, WIDTH LESS THAN 22 INCHES, ANY HEIGHT, INCLUDING ANY TYPE MOUNTING HARDWARE
- **Deep (4") Contour**  
**E2615** - POSITIONING WHEELCHAIR BACK CUSHION, POSTERIOR-LATERAL, WIDTH LESS THAN 22 INCHES, ANY HEIGHT, INCLUDING ANY TYPE MOUNTING HARDWARE
- **Extra Deep (6") Contour**  
**E2620** - POSITIONING WHEELCHAIR BACK CUSHION, PLANAR BACK WITH LATERAL SUPPORTS, WIDTH LESS THAN 22 INCHES, ANY HEIGHT, INCLUDING ANY TYPE MOUNTING HARDWARE



[stealthproducts.com](http://stealthproducts.com)

[sales@quantumrehab.co.uk](mailto:sales@quantumrehab.co.uk)

or by phone:

01869 324600