

**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk



**Dealer February 2021**



**Edge 3 Stretto with Contoured Seat**

Dealer Name: \_\_\_\_\_ Order Date: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 Phone \_\_\_\_\_ Email: \_\_\_\_\_  
 PO No: \_\_\_\_\_ Order Ref: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 Post Code: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 Post Code: \_\_\_\_\_  
 Quantum Specialist: \_\_\_\_\_

**Standard Features**

- 17 Colour Options
- Max User Weight: 21 st 6lbs. (136 kg)
- Contoured Seat with Folding Back
- Flip-up Armrests
- 2" x 5" Front Castors, Flat Free
- 2" x 6" Rear Castors, Flat Free
- SRS Suspension
- 2" x 12.5" Flat Free Drive Tyres
- NE+ Control System
- 70° Swing-Away Legrests
- Crash Tested
- 4 Mph (6 Km/h) Motors

**Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature**

**Base Model**

**Edge 3 Stretto with Contoured Seat - £3,495.00**

**Drive Wheels**

- |  |          |
|--|----------|
| <input type="checkbox"/> 2" x 12.5" Flat Free, Black Tyres | Standard |
| <input type="checkbox"/> 2" x 12.5" Pneumatic, Black Tyres | FOC      |
| <input type="checkbox"/> 3" x 14" Flat Free, Black Tyres   | £50.00   |
| <input type="checkbox"/> 3" x 14" Pneumatic, Black Tyres   | £50.00   |

**Contoured Seat Size Packages**

**Includes headrest and pelvic seat belt**

- 16"W x 16"D                       18"W x 18"D                       20"W x 20"D

**Colours**



*\* This colour choice may increase lead time – Please check with Q-rep before ordering*

## Seat Cushions

**Size provided will match seat size package**

- |                          |  |          |
|--------------------------|--|----------|
| <input type="checkbox"/> | Contoured Seat Cushion: <i>General use</i>   | Standard |
| <input type="checkbox"/> | Stealth Simplicity G: <i>General use, Fluid resistant Stealth-Tek cover - 25 stone weight limit</i>                                  | £60.00   |
| <input type="checkbox"/> | Stealth TRU-Comfort 2 SPP: <i>Skin protection, positioning, 4-way stretch Coolcore cover - 32 stone weight limit</i>                 | £140.00  |
| <input type="checkbox"/> | Stealth Solution® SPP: <i>Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 25 stone weight limit</i> | £170.00  |
| <input type="checkbox"/> | Stealth Zen SP: <i>Skin Protection, 4-way stretch, Coolcore cover - 79.3 to 181.4 kg weight limit</i>                                | £130.00  |
| <input type="checkbox"/> | Omit Seat Cushion  |          |

## Seat Positioning Adjustments

**Seat to Floor Height** - Measured from front edge of seat pan

- |  |   |   |
|--|---|---|
| <input type="checkbox"/> 16.5" Flat seat pan | <input type="checkbox"/> 17.5" with 5° posterior tilt | <input type="checkbox"/> 18.25" with 10° posterior tilt |
|--|---|---|

## Armrests

- |                          |   |          |
|--------------------------|---|----------|
| <input type="checkbox"/> | Flip-Up Height Adjustable Removable <i>Range of adjustment 11" to 14.5" from seat pan</i> | Standard |
| <input type="checkbox"/> | TB Flex Drop-In T-Arm <i>Range of adjustment 11" to 14.5" from seat pan</i>               | £85.00   |
| <input type="checkbox"/> | Single Post Height Adjustable <i>Range of adjustment 11" to 17" from seat pan</i>         | £110.00  |
| <input type="checkbox"/> | 10" Standard  |          |
| <input type="checkbox"/> | 14" Standard  |          |
| <input type="checkbox"/> | 10" Waterfall<br><small>Only available with single post height adjustable</small>         |          |
| <input type="checkbox"/> | 14" Waterfall<br><small>Only available with single post height adjustable</small>         |          |

## Non-Expandable Controls

- |                          |               |          |
|--------------------------|---------------|----------|
| <input type="checkbox"/> | NE + Controls | Standard |
| <input type="checkbox"/> | VR2 Controls  | FOC      |

## Hand Control Mounting Options

- |                          |             |          |
|--------------------------|-------------|----------|
| <input type="checkbox"/> | Fixed       | Standard |
| <input type="checkbox"/> | Swing Away  | £75.00   |
| <input type="checkbox"/> | Left Mount  |          |
| <input type="checkbox"/> | Right Mount |          |

## Legrests

**Enter seat pan to heel measurement [.....]. If no measurement is supplied, the standard factory option will be fitted**

- |                          |  |          |
|--------------------------|--|----------|
| <input type="checkbox"/> | 70° Swing-away Legrests <i>Range of adjustment 9" to 15" from seat pan to heel</i>   | Standard |
| <input type="checkbox"/> | 60° to 90° Adjustable Swing-away Legrests <i>Range of adjustment 11" to 16.5" from seat pan to heel</i>                        | £350.00  |
| <input type="checkbox"/> | Manual Elevating Swing-Away Legrests <i>With calf pads. Range of adjustment 9" to 20" from seat pan to heel</i>                | £130.00  |
| <input type="checkbox"/> | Centre Mount Foot Platform <i>With calf pads. Tool adjustable 80° to 90° . Range of adjustment 10" to 16" seat pan to heel</i> | £350.00  |

## Speed Program and Lights Package

- |                          |   |          |
|--------------------------|---|----------|
| <input type="checkbox"/> | 4 Mph Program <i>Lights and indicators not included</i> | Standard |
| <input type="checkbox"/> | 4 Mph Program with Lights and Indicators                | £450.00  |
| <input type="checkbox"/> | 6 Mph Program with Lights and Indicators                | £450.00  |

## Accessories

- |                          |  |         |
|--------------------------|--|---------|
| <input type="checkbox"/> | Fender Lights Includes seperate switch for operation. Does not include indicators. Not road legal. | £195.00 |
| <input type="checkbox"/> | Swing-Away Tray Table  | £350.00 |

## Price

### Notes

1. All measurements given are approximate.
2. By faxing or mailing this order you accept all responsibility for the configuration ordered.
3. Errors made in filling out this order, may result in further charges or delays.
4. If this order has a deadline you must advise us ASAP. We will try to deliver for this date, but cannot accept liability if we are unable to fulfil this order by your deadline.
5. Information relating to crash testing is correct at time of printing.  
Pride Mobility constantly reviews and updates this information, if you require further information about crash testing, please contact your Quantum Specialist

<b>Sub Total:</b>	
<b>Options:</b>	
<b>Total:</b>	