

**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk



**QUANTUM  
 KOZMO  
 PORTABLE**



**QUANTUM®**

**Dealer April 2021**



**Kozmo Portable**

Dealer Name: \_\_\_\_\_ Order Date: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 Post Code: \_\_\_\_\_  
 Phone \_\_\_\_\_ Email: \_\_\_\_\_  
 PO No: \_\_\_\_\_ Order Ref: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_  
 \_\_\_\_\_  
 Post Code: \_\_\_\_\_  
 Quantum Specialist: \_\_\_\_\_

**Standard Features**

- Available Seat Widths = 12" - 18"
- Available Seat Depths = 12" - 18"
- 2 x 18Ah batteries
- GC3 Joystick
- Max User Weight 17 Stone 8lbs (113 kg)
- 6" Flat Free Front Castors
- 9" Flat Free Drive Tyres
- 25.4" Turning Radius
- 3.6 Mph (5.79 Km/h) Motors
- Dismantles into 4 pieces
- Quick Height Adjustable Armrests
- Flip-up foot Platform
- 6 Base Colours
- Padded Backrest
- Auto Buckle Lap Belt, Push Button Release

**Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature**

**Base Model**

**Kozmo - £2,360.00**

**Drive Wheels**

- 9" Flat Free, Black Tyres Standard

**Castor Wheels**

- 6" Flat Free, Black Tyres Standard

**Colours**

-  Ruby Red    
   Pearl White    
   Rose Quartz    
   Citrine Yellow    
   Amber Orange    
   Sapphire Blue

**Seat Dimensions**

Select width and depth package and enter desired seat width and depth in that package range. Available in 1" increments.

**Seat Width Packages**

- 12"W-16"W      16"W-18"W  
 Set Width at: \_\_\_\_\_

**Seat Depth Packages**

- 12"D-16"D      16"D-18"D  
 Set Depth at: \_\_\_\_\_

**Seat Heights**

- Standard Seat Post For 19" or 20" seat pan to floor heights  
 Short Seat Post For 16" or 17" seat pan to floor heights  
**Seat to Floor Heights**    
  16"      17"      19"      20"

Kozmo - DEALER

