

# Quantum®

32 Wedgwood Road, Bicester,  
Oxon, OX26 4UL  
Phone: 01869 324600  
Fax: 01869 323070  
Email: sales@quantumrehab.co.uk  
Web: www.quantumrehab.co.uk



QUANTUM®

QUANTUM  
**KOZMO**  
PORTABLE



Dealer December 2018 (V1.4)



## Kozmo Portable

Dealer Name: \_\_\_\_\_ Order Date: \_\_\_\_\_

Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_ Post Code: \_\_\_\_\_

Phone \_\_\_\_\_ Email: \_\_\_\_\_

PO No: \_\_\_\_\_ Order Ref: \_\_\_\_\_

Delivery Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_ Post Code: \_\_\_\_\_

Quantum Specialist: \_\_\_\_\_

## Standard Features

- Available Seat Widths = 12" - 18"
- Available Seat Depths = 12" - 18"
- 2 x 18Ah batteries
- GC3 Joystick
- Max User Weight 17 Stone 8lbs (113 kg)
- 6" Flat Free Front Castors
- 9" Flat Free Drive Tyres
- 25.4" Turning Radius
- 3.6 Mph (5.79 Km/h) Motors
- Dismantles into 4 pieces
- Quick Height Adjustable Armrests
- Flip-up foot Platform
- 6 Base Colours
- Padded Backrest
- Auto Buckle Lap Belt, Push Button Release

Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature

## Base Model

Kozmo - £2,250.00

## Drive Wheels

- 9" Flat Free, Black Tyres Standard

## Castor Wheels

- 6" Flat Free, Black Tyres Standard

## Colours



Ruby Red



Pearl White



Rose Quartz



Citrine Yellow



Amber Orange



Sapphire Blue

## Seat Dimensions

Select width and depth package and enter desired seat width and depth in that package range. Available in 1" increments.

### Seat Width Packages

- 12"W-16"W       16"W-18"W

Set Width at: \_\_\_\_\_

### Seat Depth Packages

- 12"D-16"D       16"D-18"D

Set Depth at: \_\_\_\_\_

## Seat Heights

Standard Seat Post For 19" or 20" seat pan to floor heights

Short Seat Post For 16" or 17" seat pan to floor heights

### Seat to Floor Heights

- 16"       17"       19"       20"

## Backrest Upholstery

Visco padded folding backrest Standard

**Backrest Heights** *Seat pan to top of back*

16"     17"     18"     19"     20"

## Seat Cushions

No Seat Cushion Required FOC

Stealth Solution® SPP: *Skin protection, positioning, 4-way stretch, Coolcore cover and urethane liner - 25 stone weight limit* Standard

**Please refer to our Stealth Cushion Order form for available sizes and enter your preferred size here.....W x .....D**

## Armrests

Single Post, Removable, Quick Height Adj. Armrests, Pair *Range of adjustment 11" - 17" in 1/2" increments* Standard

Flip-up Armrests *Mounted to Back Canes. Range of adjustment 2"* FOC

Paediatric, Removable, Height Adj. Armrests *Range of adjustment 7.75" to 9.75" in .1/2" increments* FOC

**Armpads** *14" Full length armpads not recommended with seat depths of 14" or less*

Straight, 10" Desk Length Armpad, Pair  Waterfall, 10" Desk Length Armpad, Pair

Straight, 14" Full Length Armpad, Pair  Waterfall, 14" Full Length Armpad, Pair

## Headrests

Quantum Headrest *10"W x 5"H* £140.00

Stealth Small Contoured *9.75"W x 3"H* £200.00

Stealth 10" Comfort Plus *10"W x 6"H* £200.00

## Legrests

Flip-up Foot platform Standard

70° Swing-away Legrests with Flip-up Footplates *Not available with a seat width less than 16"* £110.00

Swing-away Manual Elevating Legrests with Flip-up Footplates *Not available with a seat width less than 16"* £200.00

Paediatric Fixed Foot platform £150.00

90° Paediatric Swing-Away Legrests *Range of adjustment 7"-10"L. Includes 4" W x 6"D angle & depth adj. footplates with heel loops* £310.00

## Joystick Options

GC3 Joystick Standard

Fixed, Inline Standard

Swing-Away £75.00

### Mounting Position

Left                      or                       Right

## Lateral Supports

**Stealth Thoracic Lateral Supports** *Swing-Away laterals may interfere with flip-up armrests*

Thin Hardware, Short Link, Pair £350.00

Thin Hardware, Standard Link, Pair £350.00

### Lateral Pads

3" H x 5" L Curved Lateral Pad, Pair £70.00

4" H x 6" L Curved Lateral Pad, Pair £70.00

## Price

### Notes

- All measurements given are approximate.
- By faxing or mailing this order you accept all responsibility for the configuration ordered.
- Errors made in filling out this order, may result in further charges or delays.
- If this order has a deadline you must advise us ASAP. We will try to deliver for this date, but cannot accept liability if we are unable to fulfil this order by your deadline.
- Information relating to crash testing is correct at time of printing. Pride Mobility constantly reviews and updates this information, if you require further information about crash testing, please contact your Quantum Specialist

<b>Sub Total:</b>	£2,250.00
<b>Options:</b>	
<b>Total:</b>	

## Special Instructions