

Quantum®

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk



Q1450



Q1450 with Synergy® Seating

Dealer April 2021



Q1450 with VR2 Controls

Dealer Name: _____ Order Date: _____
 Address: _____

 Phone _____ Email: _____
 PO No: _____ Order Ref: _____
 Delivery Address: _____

 Post Code: _____
 Quantum Specialist: _____

Standard Features

- Contoured Seating System
- Overall Length 38.25"
- Base Width 29.25"
- 2 x 70ah batteries
- VR2 Controls
- 3" x 9" Front Castors, Flat Free
- 4" x 14" Solid Drive Tyres
- 4 Colour Options
- 4 Mph (6 Km/h) Motors
- High-performance Motor Package
- 26.5" Turning Radius
- Heavy Duty, Base Mounted Foot platform
- Max User Weight 42 Stone 9lbs (272 kg)
- 4" Anti-tip Wheels
- 8 Amp Off Board Charger

Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature

Base Model

Q1450 with Contoured Seating and VR2 Controls - £8,600.00

Drive Wheels

3" x 14" Flat Free, Black Tyres Standard

Castor Wheels

3" x 9" Flat Free, Rear Black Tyres Standard

Colours

Standard

Special Order



Silver



Black



Red



Blue

Contoured Seating

Pinchless Hinge High Back Seat:

- Black Vinyl, High Back W/Headrest, Pinchless Hinges and Contoured Seat Cushion Standard
- 24"W x 20" D 24"W x 22" D 24"W x 24" D 24"W x 26" D
- Solid Seat Pan Only *No Seat Cushion Supplied* FOC

Contoured Seating (Continued)

Heavy Duty High Back Seat: Includes a Lumber Cut out, Static Hinges and an HD Synergy Seat Frame

Black Vinyl, High Back W/Headrest, Static Hinges and Contoured Seat Cushion £770.00

Standard Sizes

26"W x 20" D 28"W x 20" D 30"W x 22" D 32"W x 22" D

Special Sizes *May increase delivery time*

26"W x 22" D 26"W x 24" D 28"W x 22" D 28"W x 24" D
 30"W x 20" D 30"W x 24" D 32"W x 20" D 32"W x 24" D

Solid Seat Pan Only *No Seat Cushion Supplied* FOC

Contoured Seating Armrests

Armrest Type

Flip-up, Height Adjustable, Removable *Only available with Pinchless High Back Seat. Not available with HD High Back Seat* £130.00

Heavy Duty, Bariatric, Height Adjustable, Removable *Bariatric Armrests are Full Length and require Bariatric Armpads, Only available on HD High Back Seating* Standard

8"-12"H 10"-14"H

Single Post, Height Adjustable Removable *Height Adjustable from 11"-17" in 1/2" increments. Only available with HD High Back Seating* £130.00

Armpads (Pair)

Full Length Desk Length *Not available with Bariatric, Height Adjustable, Removable Armrests*

Straight Waterfall

Black Molded, Bariatric Armpads £154.00

Synergy® Seating System

Heavy Duty Synergy Seating System *Includes Synergy® Visco Backrest and Synergy® Spectrum Foam Seat Cushion* £1,200.00

Synergy® Seating Power Option

TRU-Balance® 2 HD Tilt 40° Centre of Gravity Power Tilt. 42 Stone 9lbs Weight Capacity. *Not available with Contoured Seating* £1,600.00

Synergy® Seat Dimensions

Select width and depth package and enter desired seat width and depth in that package range, available in 1" increments. Maximum seat size is 32"W x 28"D.

Seat Width Packages

20"W-24"W 24"W-26"W 26"W-28"W

28"W-30"W 30"W-32"W

Set Width at: _____ (Max seat width 32")

Seat Depth Packages

15"D-20"D 20"D-24"D 24"D-28"D

Set Depth at: _____ (Max seat depth 28")

Synergy® Back Cane Options

Synergy® Static Backrest *Tool angle adjustable from 85° to 120° in 5° increments* Standard

Back Cane Options:

1" Endo *Reduces backrest width by 2"* £40.00

1" Meso *Increases backrest width by 2"* £40.00

Back Cane Height:

16"H 17"H 20"H 21"H 24"H 25"H

Synergy® Backrest Options

Synergy Back *Backrest Height will be the same as the Back Canes Selected Above. Overall Height = Seat Pan to Top of Backrest*

Synergy Back with 1.5" Thick Back Cushion Standard

Sling Back *Backrest Height will be the same as the Back Canes Selected Above. Overall Height = Seat Pan to Top of Backrest*

Sling Back 23" Width or less £60.00

HD Sling Back 24" Width or Greater £120.00

Omit Back *Supplied with back canes only* FOC

Synergy® Armrests

Armrest Type

Bariatric, Removable, Height Adjustable (Pair) *Bariatric Armrests are Full Length and require Bariatric Armpads* Standard

8"-12"H 10"-14"H

Removable, Single Post, Height Adjustable Heavy Duty (Pair) £120.00

Height Adjustable Cantilever (Pair) *Height Adjustable from 11"-17" depending on cane height* £150.00

Synergy® Armrest Pad Options

<input type="checkbox"/> 14"L x 2"W Straight, Pair	<input type="checkbox"/> 10"L x 2"W Straight, Pair	Standard
<input type="checkbox"/> 14"L x 3"W Waterfall, Pair	<input type="checkbox"/> 10"L x 3"W Waterfall, Pair	FOC
<input type="checkbox"/> Flat, Bariatric <i>Bariatric Armrests are full length and require Bariatric Armrests</i>		FOC
<input type="checkbox"/> 14"L x 3.5"W Gel Armpad, Pair		£160.00
<input type="checkbox"/> 10"L x 3.5"W Gel Armpad, Pair		£130.00

Synergy® Headrest Options

<input type="checkbox"/> Quantum 10" W Headrest	£120.00
<input type="checkbox"/> Stealth Small Contoured 9.75"W x 3"H	£150.00
<input type="checkbox"/> Stealth Large Contoured 11"W x 6"H	£150.00
<input type="checkbox"/> Stealth 10" Comfort Plus	£180.00
<input type="checkbox"/> Stealth 14" Comfort Plus	£180.00

Please visit www.stealthproducts.com for more headrest options

Centre Mount Foot Platforms

<input type="checkbox"/> Base Mounted Foot Platform <i>Mounted to chassis, will not move with Power Tilt</i>	Standard
<input type="checkbox"/> Adjustable Centre Mount Foot Platform <i>Mounted to Seat</i>	£244.00

Footplates

<input type="checkbox"/> 10"W x 11"D Split Footplates	or	<input type="checkbox"/> 17"W x 13"D One piece Footplate	or	<input type="checkbox"/> 17"W x 16"D One piece Footplate	Standard
---	----	--	----	--	----------

Lower Extensions *Seat pan to Heel*

<input type="checkbox"/> 10-16"L <i>Seat pan to heel</i>	or	<input type="checkbox"/> 13-19"L <i>Seat pan to heel</i>	Standard
--	----	--	----------

Calf Pads


<input type="checkbox"/> Small 6"H x 6"W	or	<input type="checkbox"/> Large 8"H x 8"W, tapered to 6"W	£40.00
--	----	--	--------

Calf Panel

<input type="checkbox"/> Small 8"H x 18"W	or	<input type="checkbox"/> Large 8"H x 24"W <i>Only available on 20" seat width or greater.</i>	£84.00
---	----	---	--------

<input type="checkbox"/> Calf Panel Extension (each)	£30.00
--	--------


Footplate Extension Kits



<input type="checkbox"/> 4" Width Extension Kit	£16.00
---	--------



<input type="checkbox"/> 2" Bell Shaped, Length Extension Kit	£16.00
---	--------



<input type="checkbox"/> 6" Side Guard Kit	£28.00
--	--------

Accessories:

<input type="checkbox"/> Wheel Bumpers, Pair	£5.00
--	-------

<input type="checkbox"/> Footplate Connector <i>Not required with one piece footplate</i>	£10.00
---	--------

Individual Legrests

Choose instead of centre mount foot platform options

<input type="checkbox"/> Heavy duty, drop in, swing-away legrests with angle adjustable footplates <i>Includes heel loops</i>	£340.00
---	---------

Legrest Angle

Footplate Selection

<input type="checkbox"/> 70°	<input type="checkbox"/> 80°	<input type="checkbox"/> 6"W x 8"D	<input type="checkbox"/> 8"W x 11"D <i>Other sizes available</i>
------------------------------	------------------------------	------------------------------------	--

Please provide the seat pan to heel measurement required, minimum 12", maximum 21": _____

Joystick Handles

- | | | |
|---|---|--|
| <input type="checkbox"/> Standard Carrot Handle | | Standard |
| <input type="checkbox"/> Mushroom Handle | <input type="checkbox"/> Extended I-Handle | £35.00 (each) |
| <input type="checkbox"/> T-Bar | <input type="checkbox"/> Soft Ball Handle | <input type="checkbox"/> Bodypoint 3 Inch Goal Post Handle |
| <input type="checkbox"/> Chin Cup | <input type="checkbox"/> Medium Ball Handle | <input type="checkbox"/> Bodypoint 4 Inch Goal Post Handle |

Joystick Mounting Brackets

- | | |
|--|----------|
| <input type="checkbox"/> Fixed, Inline | Standard |
| <input type="checkbox"/> Swing-Away | £75.00 |
| <input type="checkbox"/> Multi-Axis | £125.00 |
| <input type="checkbox"/> Retract 4 Swing-Away Joystick Arm Kit <i>Full 180° Movement</i> | £150.00 |

Mounting Position

- | | | |
|---------------------------------|----|--|
| <input type="checkbox"/> Left | or | <input type="checkbox"/> Right |
| <input type="checkbox"/> Inline | or | <input type="checkbox"/> Offset <i>Not available with Fixed Joystick Bracket</i> |

Speed Program and Lights Package

- | | |
|---|----------|
| <input checked="" type="checkbox"/> 4 Mph program | Standard |
| <input type="checkbox"/> LED Light Package | POA |

Quantum Lapbelts

- | | | | |
|---------------------------------------|---------------------------------------|---------------------------------------|----------|
| <input type="checkbox"/> 50" Lap Belt | <input type="checkbox"/> 70" Lap Belt | <input type="checkbox"/> 90" Lap Belt | Standard |
| <input type="checkbox"/> 60" Lap Belt | <input type="checkbox"/> 80" Lap Belt | | |

Seat to Floor Height Options

Seat to Floor Height

- | | | | |
|---|--|--|--|
| Contoured Seating <i>inc. Seat Cushion</i> | <input type="checkbox"/> 20.5" Seat to Floor | <input type="checkbox"/> 21.5" Seat to Floor | <input type="checkbox"/> 22.5" Seat to Floor |
| Contoured Seating <i>With Seat Pan only</i> | <input type="checkbox"/> 17.5" Seat to Floor | <input type="checkbox"/> 18.5" Seat to Floor | <input type="checkbox"/> 19.5" Seat to Floor |
| HD Contoured Seating <i>26" W and wider</i> | <input type="checkbox"/> 21.5" Seat to Floor | <input type="checkbox"/> 22.5" Seat to Floor | <input type="checkbox"/> 23.5" Seat to Floor |
| HD Contoured Seating <i>With Seat Pan only, 26" W and Wider</i> | <input type="checkbox"/> 18.5" Seat to Floor | <input type="checkbox"/> 19.5" Seat to Floor | <input type="checkbox"/> 20.5" Seat to Floor |
| Synergy® Static Seat <i>From Seat Pan</i> | <input type="checkbox"/> 16.5" Seat to Floor | <input type="checkbox"/> 17.5" Seat to Floor | <input type="checkbox"/> 18.5" Seat to Floor |
| Synergy® Seating W/ HD Tilt <i>From Seat Pan</i> | <input type="checkbox"/> 18.125" Seat to Floor | <input type="checkbox"/> 19.125" Seat to Floor | <input type="checkbox"/> 20.125" Seat to Floor |

Price

Sub Total:	£8,600.00
Options:	
Total:	

Special Instructions

For additional seating options such as thoracic supports and alternative drive controls please refer to secondary order forms, entitled 'Additional Seating Options' or 'Specialty Control Input Devices'.

Quantum Notes:

1. This product not suitable for occupied transportation.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3 %. Overall measurements will vary based on seating and accessory selections.