

**Quantum®**  
 32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk



**Dealer April 2021**



**Q4 with TB Flex Seating**

Dealer Name: \_\_\_\_\_ Order Date: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 Post Code: \_\_\_\_\_  
 Phone \_\_\_\_\_ Email: \_\_\_\_\_  
 PO No: \_\_\_\_\_ Order Ref: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_  
 \_\_\_\_\_  
 Post Code: \_\_\_\_\_  
 Quantum Specialist: \_\_\_\_\_

**Standard Features**

- TB Flex Seating System
- Max User Weight: 21 st 6lbs. (136 kg)
- 2 x 50ah batteries
- NE+ Control System
- 2" x 6" Front & Rear Castors, Flat Free
- 14" Flat Free Drive Tyres
- ATX Suspension
- 4 Mph (6 Km/h) Motors
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Flip Forward Backrest
- Tie-Down System for Transportation

**Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature**

**Base Model**

**Q4 with TB Flex Seating - £3,750.00**

**Drive Wheels**

- 3" x 14" Flat Free, Black Tyres Standard
- 3" x 14" Pneumatic, Black Tyres FOC

**Castor Wheels**

- 2" x 6" Flat Free, Black Tyres, Front and Rear FOC

**Colours**

-   
 Red Corvette
-   
 Midnight Blue
-   
 Back in Black

**Seat Dimensions**

- |                              |                              |   |                              |                              |                              |                              |                              |
|------------------------------|------------------------------|---|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| <b>Seat Width</b>            |                              |   |                              | <b>Seat Depth</b>            |                              |                              |                              |
| <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17"  | <input type="checkbox"/> 18" | <input type="checkbox"/> 15" | <input type="checkbox"/> 16" | <input type="checkbox"/> 17" | <input type="checkbox"/> 18" |
| <input type="checkbox"/> 19" | <input type="checkbox"/> 20" | <input type="checkbox"/> 21" <small>Only available with tension adjustable backrest</small> |                              | <input type="checkbox"/> 19" | <input type="checkbox"/> 20" |                              |                              |

*Please note the smallest seat cushion available is 16" x 16", and the largest is 20" x 20". Please refer to our Stealth Cushion Order Form for alternative models and sizes.*

**Seat Cushions**

- Sitrite Premier Standard
- Sitrite Premier Plus £95.00
- Tru-Comfort 2 SPP £95.00

**Cushion Size**

- 16"W x 16"D
- 16"W x 18"D
- 16"W x 20"D
- 18"W x 16"D
- 18"W x 18"D
- 18"W x 20"D
- 20"W x 16"D
- 20"W x 18"D
- 20"W x 20"D

Q4 TB FLEX SEATING

## Seat Positioning Adjustments

- |  |          |
|--|----------|
| <input type="checkbox"/> Manual Tilt 0° to 30° <i>Tool adjustable with a telescopic tube</i>   | Standard |
| <input type="checkbox"/> TB-Flex® Power Tilt 0° to 30° <i>Full Drive Lock at 20° when combined with Power Recline. Uses 1 power option</i> | £300.00  |

### Seat to Floor Height

- |                              |  |   |
|------------------------------|--|---|
| <input type="checkbox"/> 16" | <input type="checkbox"/> 17" <i>Changes Tilt range to -5° to 25°</i> | <input type="checkbox"/> 18" <i>Changes Tilt range to -10° to 20°</i> |
|------------------------------|--|---|

## Backrest Upholstery



- Tension Adjustable Back**  
FOC



- Sitrite Premier**  
Standard



- Sitrite Premier Plus**  
£170.00



- ADI**  
£.....  
Part No:  
.....

**Range of Adjustment:** *Sitrite Premier and Premier Plus = 24.5" High to 28" High. (Factory Set at 24.5") Tension Adjustable = 19" to 23"*

## Backrest Positioning Adjustments

- |   |          |
|---|----------|
| <input type="checkbox"/> Manual Settable / Flip Forward <i>Range of adjustments 95° to 125°</i>   | Standard |
| <input type="checkbox"/> Attendant Controlled Manual Recline <i>Range of adjustments 95° to 125°. Not suitable for occupied transportation.</i> | £150.00  |
| <input type="checkbox"/> Powered Recline <i>Range of adjustments 95° to 125°. Uses 1 power option.</i>  | £200.00  |

## Armrests

- |   |          |
|---|----------|
| <input type="checkbox"/> TB Flex Removable Armrest with Quick Release <i>Range of adjustments 9" to 14"</i> | Standard |
| <input type="checkbox"/> TB Flex Removable Armrest with Armrest Lever <i>Range of adjustments 9" to 14"</i> | £85.00   |
| <input type="checkbox"/> 2-Post Flip Up Removable Armrest <i>Range of adjustments 9" to 13"</i>             | £280.00  |

### Armpads

- |                                       |                                       |   |   |
|---------------------------------------|---------------------------------------|---|---|
| <input type="checkbox"/> 10" Standard | <input type="checkbox"/> 14" Standard | <input type="checkbox"/> 10" Waterfall<br><i>Only available with 2-post flip up</i> | <input type="checkbox"/> 14" Waterfall<br><i>Only available with 2-post flip up</i> |
|---------------------------------------|---------------------------------------|---|---|

## Non-Expandable Controls

- |  |          |
|--|----------|
| <input type="checkbox"/> NE + Controls <i>2 power option ready</i> | Standard |
| <input type="checkbox"/> VR2 Controls <i>2 power option ready</i>  | FOC      |

## Joystick Mounting Options

- |   |          |  |         |
|---|----------|--|---------|
| <input type="checkbox"/> Fixed                  | Standard | <input type="checkbox"/> Retract 4 Swing-away Joystick Arm Kit <i>(Full 180° Movement)</i> | £150.00 |
| <input type="checkbox"/> Swing-Away             | £75.00   | <input type="checkbox"/> Dual Control  | £380.00 |
| <input type="checkbox"/> Attendant Control Only | £100.00  |  |         |

### Mounting Position

- Left      or       Right

## Legrests

- |   |          |
|---|----------|
| <input type="checkbox"/> 70° Swing-away Legrests <i>Range of adjustments 12" to 15" from seat pan to heel</i>   | Standard |
| <input type="checkbox"/> Base Mounted Foot Platform <i>Mounted to Chassis, will not move with Power Tilt</i>  | £75.00   |
| <input type="checkbox"/> Synergy Centre Mount Foot Platform <i>Mounted to Seat. Range of adjustments 10" to 16" from seat pan to heel</i>                                     | £315.00  |
| <input type="checkbox"/> Quantum 60° to 90° Adjustable Swing-away Legrests <input type="checkbox"/> 11"-13" <input type="checkbox"/> 13.5"-16.5" <i>from seat pan to heel</i> | £337.00  |
| <input type="checkbox"/> Manual Elevating Swing-Away Legrests <i>Range of adjustments 15" to 20" from seat pan to heel</i>  | £250.00  |
| <input type="checkbox"/> Powered Elevating Swing-Away Legrests <i>Range of adjustments 15" to 20" from seat pan to heel. Uses 2 power options</i>                             | £650.00  |

## Headrests

- |   |         |
|---|---------|
| <input type="checkbox"/> Stealth 10" Comfort Plus CP250             | £140.00 |
| <input type="checkbox"/> Stealth 14" Comfort Plus CP150             | £140.00 |
| <input type="checkbox"/> Stealth 14" Adjustable Comfort Plus ACP150 | £250.00 |
| <input type="checkbox"/> Stealth Small Contoured STL-OTSM250        | £200.00 |
| <input type="checkbox"/> Stealth Large Contoured STL-OTLG150        | £200.00 |

*Please visit [www.stealthproducts.com](http://www.stealthproducts.com) for more headrest options*

## Lateral and Trunk Supports

### Stealth Thoracic Lateral Supports

#### Left Lateral Package

Lateral Thoracic Link Hardware , Cane Mount & Pad      £195.00

#### Left Lateral Pads Choose one option only

3" H x 5" L Curved Lateral Pad

4" H x 6" L Curved Lateral Pad

#### Right Lateral Package

Lateral Thoracic Link Hardware , Cane Mount & Pad      £195.00

#### Right Lateral Pads Choose one option only

3" H x 5" L Curved Lateral Pad

4" H x 6" L Curved Lateral Pad

## Pelvic Supports

Auto Buckle Lap Belt, Push Button Release Standard

#### i-Fit Pelvic Belts Supplied un-installed

Padded Offers a large padded surface area £70.00

Non-Padded Narrowed section for smaller support surfaces. Not available with 4-point £65.00

#### Size

Extra Small  
Not available with non-padded

Small

Medium

Large

Extra Large  
Not available with non-padded or 4-point.

#### Attachment Method

2-Point

4-Point Single pull only

#### Release Mechanism

Push Button

Side Release

#### Adjustment

Single Pull

Dual Pull

## Accessories

Quick Release Transfer Bars Not Compatible with 2-Post Flip Up Removable Armrest £125.00

Quantum Backpack £40.00

Swing-Away Tray Table £350.00

## Price

|                   |           |
|-------------------|-----------|
| <b>Sub Total:</b> | £3,750.00 |
| <b>Options:</b>   |           |
| <b>Total:</b>     |           |

### Notes

1. All measurements given are approximate.
2. By faxing or mailing this order you accept all responsibility for the configuration ordered.
3. Errors made in filling out this order, may result in further charges or delays.
4. If this order has a deadline you must advise us ASAP. We will try to deliver for this date, but cannot accept liability if we are unable to fulfil this order by your deadline.
5. Information relating to crash testing is correct at time of printing.  
Pride Mobility constantly reviews and updates this information, if you require further information about crash testing, please contact your Quantum Specialist