

Quantum®

32 Wedgwood Road, Bicester,
 Oxon, OX26 4UL
 Phone: 01869 324600
 Fax: 01869 323070
 Email: sales@quantumrehab.co.uk
 Web: www.quantumrehab.co.uk



Dealer April 2017 (v1.0)



Quantum Lightning with Synergy® Seating and R-Net Controls

Dealer Name: _____ Order Date: _____
 Address: _____

 Phone _____ Email: _____
 PO No: _____ Order Ref: _____
 Delivery Address: _____

 Post Code: _____

 Post Code: _____

 Quantum Specialist: _____

Standard Features

- Synergy Seating System
- 4" of Seat Width Adjustment
- 4" of Seat Depth Adjustment
- 2 x 70ah batteries
- R-Net Joystick
- 3" x 9" Front Castors, Pneumatic
- 3" x 14" Pneumatic Drive Tyres
- 5 Colour Options
- 4 Mph (6 Km/h) Motors
- ATX Suspension System
- Single Post Height Adj. Removable Armrests
- 60°-90° Swing-Away Legrests
- Maximum User Weight 21 Stone 6lbs
- Synergy Simplicity Seat Cushion
- Tie-Down System for Transport

Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature

Base Model

Lightning with Synergy® Seating and R-Net Controls - £3,850.00

Drive Wheels

- 3" x 14" Pneumatic, Grey Tyres Standard
- 3" x 14" Flat Free, Grey Tyres £150.00

Castor Wheels

- 3" x 9" Pneumatic, Grey Tyres Standard
- 3" x 9" Flat Free, Grey Tyres £60.00

Colours



Seat Dimensions

Select width and depth package and enter desired seat width and depth in that package range, available in 1" increments. Maximum seat size is 22"W x 22"D.

Seat Width Packages

- 14"W-18"W
 - 16"W-20"W
- Set Width at: _____ (Max seat width 20")

Seat Depth Packages

- 14"D-18"D
 - 16"D-20"D
- Set Depth at: _____ (Max seat depth 20")

Seat Cushions

- | | |
|---|----------|
| <input type="checkbox"/> No Seat Cushion Required | FOC |
| <input type="checkbox"/> Simplicity Cushion: <i>General Use</i> | Standard |
| <input type="checkbox"/> TRU-Comfort™2 4-Way Stretch: <i>Skin Protection</i> | £45.00 |
| <input type="checkbox"/> Solution® Cushion: <i>Skin Protection & Positioning</i> | £100.00 |
| <input type="checkbox"/> Solution® 1 Cushion: <i>Skin Protection</i> | £100.00 |
| <input type="checkbox"/> Spectrum® Gel Cushion: <i>Positioning</i> | £160.00 |
| <input type="checkbox"/> Spectrum® Foam Cushion: <i>Skin Protection & Positioning</i> | £160.00 |

Refer to cushion matrix for size availability and enter your selection here: W..... X D.....

Synergy® Backrest Specifications

- | | |
|---|----------|
| <input type="checkbox"/> Synergy® Backrest Assembly <i>Including Visco padded backrest</i> | Standard |
| <input type="checkbox"/> Tension Adjustable Backrest Upholstery <i>This option may increase the delivery time</i> | £120.00 |
| <input type="checkbox"/> Back Canes Only <i>Suitable for third party backrests. A backrest will not be supplied</i> | FOC |

Backrest/Back Cane Height:

Set Backrest/Back Cane Height at: _____ (Available in 1" increments. Min 13", Max 25") *Please refer to Synergy Backrest Matrix on page* for seat width and backrest height compatibility*

Back Height = Seat Pan to Top of Back | Odd heights = 4" gap at bottom, Even heights = 3" gap at bottom.

TRU-Comfort Plus® Backrest Specifications

- | | |
|---|----------|
| <input type="checkbox"/> TRU-Comfort Plus® Backrest Assembly <i>4-way stretch fabric (Positioning). Only available with seat widths 16", 18", 20" and 22"</i> | £220.00 |
| Backrest Height: | |
| <input type="checkbox"/> 27" Tall | Standard |
| <input type="checkbox"/> 24" Tall | FOC |

Back Height = Seat Pan to Top of Back | 3" gap at bottom.

Seat Positioning Adjustments

- | | |
|--|----------|
| <input type="checkbox"/> Static Seat | Standard |
| <input type="checkbox"/> TRU-Balance®2 Power Tilt <i>0-55 degree centre of gravity (CG) power tilt seating system. Uses 1 power option</i> | £1100.00 |

Seat to Floor Height

- | | | | |
|----------------|--|---|--|
| Static Seat | <input type="checkbox"/> 16.5" Seat to Floor | <input type="checkbox"/> 17.25" Seat to Floor | <input type="checkbox"/> 18" Seat to Floor |
| TB2 Power Tilt | <input type="checkbox"/> 18" Seat to Floor | <input type="checkbox"/> 18.25" Seat to Floor | <input type="checkbox"/> 19.5" Seat to Floor |

Backrest Positioning Adjustments

- | | |
|---|----------|
| <input type="checkbox"/> Static Backrest <i>Tool angle adjustable from 85° to 120° in 5° increments</i> | Standard |
| Additional Static Back Options: | |
| <input type="checkbox"/> Flip Forward Backrest <i>Not suitable for occupied transportation</i> | £200.00 |
| <input type="checkbox"/> 1" Endo Brackets <i>Reduces backrest width by 2". Not available with TRU-Comfort Plus backrest on seat widths 14"-17"</i> | £40.00 |
| <input type="checkbox"/> 2" Endo Brackets <i>Reduces backrest width by 4". Not available with TRU-Comfort Plus backrest on seat widths 14"-19"</i> | £40.00 |
| <input type="checkbox"/> TRU-Balance®2 Power Recline <i>Only available with 16"-22" s/widths. 168° of power recline with full drive lockout at 20°. Uses 1 power option. Limited backrest heights available. See reference chart on page*</i> | £1790.00 |

Armrests

- | | |
|---|----------|
| <input type="checkbox"/> Single Post, Quick Height Adjustable Removable Armrests, Pair <i>Height adjustment from 11" - 17" in 1/2" increments</i> | Standard |
| <input type="checkbox"/> Single Post, Tool Height Adjustable Armrests, Heavy Duty, Removable, Pair <i>Height adjustment from 10"-14"</i> | £220.00 |
| <input type="checkbox"/> Paediatric, Removable, Height Adjustable <i>Height adjustment from 7.75" to 9.75" in 1/2" increments. Not available with full length ampads</i> | FOC |
| <input type="checkbox"/> Cantilever Flip-Up Armrests, Pair <i>Not available with a tension adjustable backrest. Increases overall width by 7". Height adj. 11"-17" depending on backrest height. Not available with Flip Forward Backrest. Not available with TB2 Power Recline</i> | £280.00 |
| <input type="checkbox"/> 2-Post, Flip Back, Height Adjustable, Removable <i>Height range 9"-13". Not available with a seat depth less than 16"</i> | £145.00 |

Armrests Pad Options

- | | | |
|--|---|----------|
| <input type="checkbox"/> 14"L x 2"W Straight, Pair | <input type="checkbox"/> 10"L x 2"W Straight, Pair | Standard |
| <input type="checkbox"/> 14"L x 3"W Waterfall, Pair | <input type="checkbox"/> 10"L x 3"W Waterfall, Pair | FOC |
| <input type="checkbox"/> 14"L x 3.5"W Gel Armpad, Pair | | £160.00 |
| <input type="checkbox"/> 10"L x 3.5"W Gel Armpad, Pair | | £130.00 |

Ottobock Arm Channel (Each)

- | | |
|---|---------|
| <input type="checkbox"/> Medium (4.75" x 11.5") <i>Desk length armrest required</i> | £140.00 |
| <input type="checkbox"/> Large (4.75" x 14") <i>Full length armrest required</i> | £140.00 |
| <input type="checkbox"/> One Piece (4.75" x 18") <i>Full length armrest required. Not compatible with Hand Pad.</i> | £140.00 |
| <input type="checkbox"/> Ottobock Hand Pad <i>Not available on same side as controller</i> | £50.00 |
| <input type="checkbox"/> Medium <input type="checkbox"/> Large | |

Ottobock Arm Channel Mounting Position

- Left Right Both *Please amend price for pair*

Headrests

- | | |
|---|---------|
| <input type="checkbox"/> Quantum 10" W Headrest | £140.00 |
| <input type="checkbox"/> Stealth Small Contoured 9.75"W x 3"H | £200.00 |
| <input type="checkbox"/> Stealth Large Contoured 11"W x 6"H | £200.00 |
| <input type="checkbox"/> Stealth 10" Comfort Plus | £200.00 |
| <input type="checkbox"/> Stealth 14" Comfort Plus | £200.00 |

Please visit www.stealthproducts.com for more headrest options

Individual Legrests

- | | |
|---|----------|
| <input type="checkbox"/> 60-90° Swing-Away Legrests <i>Tool adjustable in 5° increments. Factory set at 70°</i> | Standard |
| Footplate Selection <i>Includes heel loops</i> | |
| <input type="checkbox"/> Flip-up, Extruded Footplates, 7.25"W x 6"D, Pair | Standard |
| <input type="checkbox"/> Angle Adjustable Footplates, 6"W x 8"D, Pair <i>Other sizes available</i> | £115.00 |

Please provide the seat pan to heel measurement required, minimum 11", maximum 19.5": _____

- | | |
|--|---|
| <input type="checkbox"/> Heavy duty, drop in, swing-away legrests with angle adjustable footplates | £130.00 |
| Legrest Angle Footplate Selection <i>Includes heel loops</i> | |
| <input type="checkbox"/> 70° <input type="checkbox"/> 80° | <input type="checkbox"/> 6"W x 8"D <input type="checkbox"/> 8"W x 11"D <i>Other sizes available</i> |

Please provide the seat pan to heel measurement required, minimum 12", maximum 18": _____

- | | |
|---|-----|
| <input type="checkbox"/> 90° Paediatric Swing-Away Legrests <i>Seat Pan to Heel range of adjustment 7"-10"L. Includes 4" x 6"D angle adjustable footplates and heel loops</i> | FOC |
|---|-----|

- | | |
|---|----------|
| <input type="checkbox"/> 70° Tapered Legrests <i>Not available with seat widths 17" or less</i> | FOC |
| Footplate Selection <i>Includes heel loops</i> | |
| <input type="checkbox"/> Flip-up, Extruded Footplates, 7.25"W x 6"D, Pair | Standard |
| <input type="checkbox"/> Angle Adjustable Footplates, 6"W x 8"D, Pair | £115.00 |

Please provide the seat pan to heel measurement required, minimum 11", maximum 19.5": _____

Elevating Legrests

Choose instead of centre mount foot platform options

- | | |
|---|----------|
| <input type="checkbox"/> Manual Elevating Swing-Away Legrests <i>Not suitable for occupied transportation</i> | £130.00 |
| Footplate Selection <i>Includes heel loops</i> | |
| <input type="checkbox"/> Flip-up, Extruded Footplates, 7.25" W x 6"D, Pair | Standard |
| <input type="checkbox"/> Angle Adjustable Footplates, 6" W x 8"D, Includes Wheel Bumpers, Pair <i>Other sizes available</i> | £115.00 |

Please provide the seat pan to heel measurement required, minimum 11", maximum 18": _____

Elevating Legrests (Continued)

- 70° Euro Power Elevating Drop-in Legrests. *Includes calf pads. Uses 2 power options.* £720.00
- Footplate Selection** *Includes heel loops*
- Flip-up, Extruded Footplates, 7.25"W x 6"D, Pair Standard
- Angle Adjustable Footplates, 6"W x 8"D, Pair *Other sizes available* £115.00

Please provide the seat pan to heel measurement required, minimum 11", maximum 18": _____

- 70° Power Articulating Elevating Swing-Away Legrests 4" articulation range. *Uses 2 power options.* £1320.00
- Footplate Selection** *Includes heel loops*
- Angle Adjustable Footplates, 6"W x 8"D, Includes Wheel Bumpers, Pair *Other sizes available* Standard

Please provide the seat pan to heel measurement required, minimum 11", maximum 18": _____

Accessories:

- Small Calf Pads, Pair 5.35"W-6"H or Large Calf Pads, Pair 5.86"Wx8"H £88.00
- Gel Padded Adductor Buttons, Pair £30.00
- Operation Method:** Independent Combined

Additional Legrest Options

- Stump Board support (each) *Includes a 60-90° Swing-Away Legrest hanger* Left Right £300.00
- Residual Limb support (each) Left Right £180.00

If both a left and right support are required, please amend price accordingly

R-Net Expandable Controls

- R-Net Joystick Standard
- Mount Joystick as Attendant Control Only £100.00
- Dual Control £430.00
- R-Net Omni Enhanced Display *Includes x2 egg switches and swing-away display mount. Does not include R-Net Joystick* £763.00

Joystick Handles

- Standard Carrot Handle Standard
- Mushroom Handle Extended I-Handle Bodypoint 3 Inch Goal Post Handle £35.00 (each)
- T-Bar Soft Ball Handle Bodypoint 4 Inch Goal Post Handle
- Chin Cup Medium Ball Handle

Joystick Mounting Brackets



- Fixed, Inline Standard
- Swing-Away £75.00
- Multi-Axis £125.00
- Retract 4 Swing-Away Joystick Arm Kit *Full 180° Movement* £150.00
- Mounting Position**
- Left or Right
- Inline or Offset *Not available with Fixed Joystick Bracket*

Speed Program and Lights Package

- 4 Mph program *Available with VR2 or R-Net. Lights and indicators not included* Standard
- 4 Mph program, with Lights and indicators £400.00

Quantum Accessories

Some accessories will be shipped with the unit, not installed. Multiple accessories can be mounted along the various tracks, some interference may occur.

-  Glove Box *Supplied uninstalled* £30.00
-  Clothing Guard (pair) *Only available with Flip Up armrests* £100.00

