

# Quantum®

32 Wedgwood Road, Bicester,  
 Oxon, OX26 4UL  
 Phone: 01869 324600  
 Fax: 01869 323070  
 Email: sales@quantumrehab.co.uk  
 Web: www.quantumrehab.co.uk



## Dealer April 2021



### Sparky with ION Seating

Dealer Name: \_\_\_\_\_ Order Date: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 Post Code: \_\_\_\_\_  
 Phone \_\_\_\_\_ Email: \_\_\_\_\_  
 PO No: \_\_\_\_\_ Order Ref: \_\_\_\_\_  
 Delivery Address: \_\_\_\_\_  
 \_\_\_\_\_  
 Post Code: \_\_\_\_\_  
 Quantum Specialist: \_\_\_\_\_

### Standard Features

- Seat Width Adjustment 10" - 14"
- Seat Depth Adjustment 10" - 14"
- 2 x 40ah batteries
- VR2 6 Key Swing-away Joystick
- Max User Weight 8 Stone 13lbs
- 8" x 2" Pneumatic Front Castors
- 14" x 4" Pneumatic Drive Tyres
- Powered Seat Tilt
- 4 Mph (6 Km/h) Motors
- Rear Suspension
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Attendant Operated Manual Recline
- Tension Adjustable Backrest
- Quantum 8" Headrest

**Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature**

### Base Model

**Sparky with ION Seating - £3,310.00**

### Drive Wheels

- 14" x 4" Pneumatic, Black Tyres Standard
- 14" x 4" Flat Free, Black Tyres £125.00

### Castor Wheels

- 8" x 2" Pneumatic, Black Tyres Standard
- 8" x 2" Flat Free, Black Tyres £55.00

### Colours



White



Blue

### Seat Dimensions

- 10" W x 10" D
- 12" W x 12" D
- 14" W x 14" D Standard

*Seat width and seat depth can be set between 10" and 14", in 1" increments e.g. 11" w x 14" d.  
 Standard seat cushion supplied as 10" x 10", 12" x 12" or 14" x 14". For alternative cushion options please contact your Quantum Specialist*

### Seat Cushions

- ION Seat Cushion with Air Mesh Supplied as 10" x 10", 12" x 12" or 14" x 14". Standard
- Synergy Solution POA

*Please contact your Quantum Specialist for alternative cushions and sizes  
 JAZZY SPARKY ION SEATING - DEALER*

## Seat Positioning Adjustments

Power Tilt 0° to 30° *Seat to ground height = 16.5" (1 Actuator)* Standard

## Backrest Upholstery



Tension Adjustable Air Mesh, 14" Tall  
Standard



Stealth ADI Backrest  
POA  
*Please contact your Quantum Specialist for ADI options*

## Backrest Positioning Adjustments

Attendant Controlled Manual Recline *Range of adjustments 95° to 125°* Standard

## Armrests

ION height adjustable, removable armrests, with 10" L pads *Range of adjustments = 2"* Standard

## Joystick Options

VR2 Joystick with Swing-away controller mount  Left  Right Standard

Mount joystick as attendant control only £100.00

Dual controls £380.00

## Legrests

70° Swing-away legrests with 4" x 6" angle and depth adjustable footplates *Includes heel loops* Standard

Swing-away manual elevating legrest £60.00 (each)

Left  Right

## Headrests

Quantum Headrest with 8" Pad Standard

Stealth Small Contoured £150.00

Stealth 10" Comfort Plus £180.00

*Please refer to Stealth catalogue, for additional options*

## Posture Supports

Auto Buckle Lap Belt, Push Button Release Standard

iFit Pelvic Belt *Please refer to iFit catalogue, for additional options* POA

iFit Pelvic Posture Support Harness *Please refer to iFit catalogue, for additional options* POA

Stealth ADI Thoracic Lateral Supports *Only compatible with ADI backrest option* POA

## Accessories

Swing-Away Tray Table £350.00

## Price

<b>Sub Total:</b>	£3,310.00
<b>Options:</b>	
<b>Total:</b>	

### Notes

1. All measurements given are approximate.
2. By faxing or mailing this order you accept all responsibility for the configuration ordered.
3. Errors made in filling out this order, may result in further charges or delays.
4. If this order has a deadline you must advise us ASAP. We will try to deliver for this date, but cannot accept liability if we are unable to fulfil this order by your deadline.
5. Information relating to crash testing is correct at time of printing.  
Pride Mobility constantly reviews and updates this information, if you require further information about crash testing, please contact your Quantum Specialist

## Special Instructions

---



---



---



---