



**QUANTUM**

**Cushion Order Form**



**Dealer December 2018 (V1.2)**

**Stealth Cushions Order Form**

Dealer Name: \_\_\_\_\_ Order Date: \_\_\_\_\_

Address: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Post Code: \_\_\_\_\_

Phone \_\_\_\_\_ Email: \_\_\_\_\_

PO No: \_\_\_\_\_ Order Ref: \_\_\_\_\_

Delivery Address: \_\_\_\_\_

\_\_\_\_\_


Post Code: \_\_\_\_\_

Quantum Specialist: \_\_\_\_\_

**Simplicity G (25 Stone)**

**Molded High Density foam base**

- for client with mild skin protection and support needs and minimal risk of skin breakdown.
- soft, gently contoured for comfort, promotes sitting tolerance
- mild trochanter support for stability



**Coccyx/Sacral relief incorporated into cushion**

- reduces surface tension in the sacral/IT well
- promotes off-loading coccyx for improved comfort and sitting tolerance

**PSIS contour /lift incorporated into rear of cushion**

- Increases surface contact area for sitting

**Fluid resistant Stealth-Tek cover**

- for hygiene, comfort
- provides improved infection control/IC barrier

**TRU-Comfort 2 SPP (32 Stone)**

**Deep contoured molded HD foam base with pelvic well for maximum positioning and skin protection needs**

- provide neutral pelvic support
- reduce anterior pelvic drift into PPT
- fatigue reduction



**1" Visco Elastic foam overlay**

- provides deeper immersion into foam for increased load distribution

**Enhanced waterfall front cushion design**

- reduces areas of potential pressure areas on posterior thighs,
- allow greater skin surface area contact and improved overall load distribution

- promotes proper positioning of lower extremities

**Mild waterfall on front of cushion**

- provides increased comfort
- reduces potential pressure points on posterior thighs

**Cool Core cover**

- enhanced foam immersion,
- moisture and temperature control for enhanced skin protection

**Solution SPP (25 Stone)**

**Contoured HD foam base**

- Superior positioning and pelvic support

**Mild lateral shelf**

- provides trochanter support for improved stability and load distribution

**Pre-ischial shelf**

- provides resistance to shearing and anterior pelvic drift

**Twin Gel insert in IT/Sacral well**

- provides superior IT immersion and offloading

**Visco-Elastic foam overlay**


- improves overall surface area distribution and pressure reduction across the entire

**Mild waterfall front edge**

- reduces potential pressure points in posterior thighs

**Cool Core cover**

- enhanced foam immersion
- moisture and temperature control for enhanced skin protection



**Spectrum Foam SPP (46 Stone)**

**Specialized pressure management and maximum comfort provided by a high density foam base with VE foam overlay**

**Visco Elastic foam overlay**

- improves overall surface area distribution and pressure reduction across the entire cushion surface



**Lateral and pre-ischial anti-thrust shelf**

- provides pelvic support and helps prevent anterior pelvic drift into posterior tilt

**Ischial VE insert in the bony area for offloading**

- allows for total immersion of ITs thus trochanters are supported on lateral support shelf to maximize stability and distribute support load

**Mild waterfall on front of cushion**

- provides increased comfort
- reduces potential pressure points on posterior thighs

**Cool Core cover**

- enhanced foam immersion
- moisture and temperature control for enhanced skin protection

**Spectrum Gel SPP Cushion (25 Stone)**

**Specialized pressure management and maximum comfort provided by a high density foam base with VE foam overlay**

**Visco Elastic foam overlay**

- improves overall surface area distribution and pressure reduction across the entire cushion surface

**Lateral and pre-ischial anti-thrust shelf**

- provides pelvic support and helps prevent anterior pelvic drift into posterior tilt

**Ischial VE insert in the bony area for offloading**




- allows for total immersion of ITs thus trochanters are supported on lateral support shelf to maximize stability and distribute support load

**Mild waterfall on front of cushion**

- provides increased comfort
- reduces potential pressure points on posterior thighs

**Cool Core cover**

- enhanced foam immersion
- moisture and temperature control for enhanced skin protection



COMING SOON

## Simplicity G Cushion

Size	Part	Price	Size	Part	Price
☐ 16"W x 16"D	ST-SIMG1616	£70.00	☐ 20"W x 16"D	ST-SIMG2016	£80.00
☐ 16"W x 18"D	ST-SIMG1618	£70.00	☐ 20"W x 18"D	ST-SIMG2018	£80.00
☐ 18"W x 16"D	ST-SIMG1816	£70.00	☐ 24"W x 18"D	ST-SIMG2016	£80.00
☐ 18"W x 18"D	ST-SIMG1818	£70.00			

## TRU-Comfort 2 SPP Cushion

Size	Part	Price	Size	Part	Price
☐ 16"W x 14"D	ST-TC2SPP1614	£150.00	☐ 20"W x 20"D	ST-TC2SPP2020	£150.00
☐ 16"W x 16"D	ST-TC2SPP1616	£150.00	☐ 20"W x 24"D	ST-TC2SPP2024	£170.00
☐ 16"W x 18"D	ST-TC2SPP1618	£150.00	☐ 22"W x 18"D	ST-TC2SPP2218	£170.00
☐ 18"W x 16"D	ST-TC2SPP1816	£150.00	☐ 22"W x 20"D	ST-TC2SPP2220	£170.00
☐ 18"W x 18"D	ST-TC2SPP1818	£150.00	☐ 22"W x 22"D	ST-TC2SPP2222	£170.00
☐ 18"W x 20"D	ST-TC2SPP1820	£150.00	☐ 22"W x 24"D	ST-TC2SPP2224	£170.00
☐ 20"W x 18"D	ST-TC2SPP2018	£150.00			

## Solution SPP Cushion (Includes Urethane Liner)

Size	Part	Price	Size	Part	Price
☐ 10"W x 10"D	ST-SOLSPP1010	£180.00	☐ 18"W x 18"D	ST-SOLSPP1818	£180.00
☐ 10"W x 12"D	ST-SOLSPP1012	£180.00	☐ 18"W x 20"D	ST-SOLSPP1820	£180.00
☐ 12"W x 10"D	ST-SOLSPP1210	£180.00	☐ 18"W x 22"D	ST-SOLSPP1822	£180.00
☐ 12"W x 12"D	ST-SOLSPP1212	£180.00	☐ 19"W x 17"D	ST-SOLSPP1917	£180.00
☐ 12"W x 14"D	ST-SOLSPP1214	£180.00	☐ 19"W x 18"D	ST-SOLSPP1918	£180.00
☐ 14"W x 12"D	ST-SOLSPP1414	£180.00	☐ 19"W x 19"D	ST-SOLSPP1919	£180.00
☐ 14"W x 14"D	ST-SOLSPP1414	£180.00	☐ 19"W x 20"D	ST-SOLSPP1920	£180.00
☐ 14"W x 16"D	ST-SOLSPP1416	£180.00	☐ 20"W x 16"D	ST-SOLSPP2016	£180.00
☐ 14"W x 18"D	ST-SOLSPP1414	£180.00	☐ 20"W x 18"D	ST-SOLSPP12018	£180.00
☐ 16"W x 16"D	ST-SOLSPP1616	£180.00	☐ 20"W x 20"D	ST-SOLSPP2020	£180.00
☐ 16"W x 17"D	ST-SOLSPP1617	£180.00	☐ 20"W x 22"D	ST-SOLSPP2022	£200.00
☐ 16"W x 18"D	ST-SOLSPP1618	£180.00	☐ 22"W x 18"D	ST-SOLSPP2218	£200.00
☐ 16"W x 20"D	ST-SOLSPP1620	£180.00	☐ 22"W x 20"D	ST-SOLSPP2220	£200.00
☐ 17"W x 17"D	ST-SOLSPP1717	£180.00	☐ 22"W x 22"D	ST-SOLSPP2222	£200.00
☐ 17"W x 19"D	ST-SOLSPP1719	£180.00	☐ 24"W x 20"D	ST-SOLSPP2420	£200.00
☐ 18"W x 16"D	ST-SOLSPP1816	£180.00	☐ 24"W x 22"D	ST-SOLSPP2422	£200.00

## Spectrum Foam SPP Cushion (Includes Urethane Liner)

Size	Part	Price	Size	Part	Price
☐ 16"W x 16"D	ST-SPCFSP1616	£250.00	☐ 19"W x 17"D	ST-SPCFSP1917	£250.00
☐ 16"W x 18"D	ST-SPCFSP1618	£250.00	☐ 19"W x 19"D	ST-SPCFSP1919	£250.00
☐ 16"W x 20"D	ST-SPCFSP1620	£250.00	☐ 20"W x 18"D	ST-SPCFSP2018	£250.00
☐ 17"W x 17"D	ST-SPCFSP1717	£250.00	☐ 20"W x 20"D	ST-SPCFSP2020	£250.00
☐ 17"W x 18"D	ST-SPCFSP1718	£250.00	☐ 20"W x 22"D	ST-SPCFSP2022	£300.00
☐ 17"W x 19"D	ST-SPCFSP1719	£250.00	☐ 20"W x 24"D	ST-SPCFSP2024	£300.00
☐ 18"W x 16"D	ST-SPCFSP1816	£250.00	☐ 22"W x 18"D	ST-SPCFSP2218	£300.00
☐ 18"W x 18"D	ST-SPCFSP1818	£250.00	☐ 22"W x 20"D	ST-SPCFSP2220	£300.00
☐ 18"W x 20"D	ST-SPCFSP1820	£250.00	☐ 22"W x 22"D	ST-SPCFSP2222	£300.00
☐ 18"W x 22"D	ST-SPCFSP1822	£250.00	☐ 22"W x 24"D	ST-SPCFSP2224	£300.00

## Spectrum Foam SPP Cushion - Continued (Includes Urethane Liner)

Size	Part	Price	Size	Part	Price
<input type="checkbox"/> 24"W x 20"D	ST-SPCFSP2420	£300.00	<input type="checkbox"/> 26"W x 32"D	ST-SPCFSP2632	£300.00
<input type="checkbox"/> 24"W x 22"D	ST-SPCFSP2422	£300.00	<input type="checkbox"/> 28"W x 24"D	ST-SPCFSP2824	£300.00
<input type="checkbox"/> 24"W x 24"D	ST-SPCFSP2424	£300.00	<input type="checkbox"/> 28"W x 26"D	ST-SPCFSP2826	£300.00
<input type="checkbox"/> 26"W x 24"D	ST-SPCFSP2624	£300.00	<input type="checkbox"/> 28"W x 28"D	ST-SPCFSP2828	£300.00
<input type="checkbox"/> 26"W x 26"D	ST-SPCFSP2626	£300.00	<input type="checkbox"/> 30"W x 28"D	ST-SPCFSP3028	£300.00
<input type="checkbox"/> 26"W x 28"D	ST-SPCFSP2628	£300.00	<input type="checkbox"/> 30"W x 30"D	ST-SPCFSP3030	£300.00
<input type="checkbox"/> 26"W x 30"D	ST-SPCFSP2630	£300.00	<input type="checkbox"/> 30"W x 32"D	ST-SPCFSP3032	£300.00
			<input type="checkbox"/> 32"W x 28"D	ST-SPCFSP3228	£300.00

## Spectrum Gel SPP Cushion (Includes Urethane Liner)

Size	Part	Price	Size	Part	Price
<input type="checkbox"/> 14"W x 18"D	ST-SPCGSP1418	£250.00	<input type="checkbox"/> 19"W x 20"D	ST-SPCGSP1920	£250.00
<input type="checkbox"/> 16"W x 16"D	ST-SPCGSP1616	£250.00	<input type="checkbox"/> 20"W x 16"D	ST-SPCGSP2016	£250.00
<input type="checkbox"/> 16"W x 18"D	ST-SPCGSP1618	£250.00	<input type="checkbox"/> 20"W x 18"D	ST-SPCGSP2018	£250.00
<input type="checkbox"/> 16"W x 19"D	ST-SPCGSP1619	£250.00	<input type="checkbox"/> 20"W x 20"D	ST-SPCGSP2020	£300.00
<input type="checkbox"/> 17"W x 17"D	ST-SPCGSP1717	£250.00	<input type="checkbox"/> 20"W x 22"D	ST-SPCGSP2022	£300.00
<input type="checkbox"/> 17"W x 18"D	ST-SPCGSP1718	£250.00	<input type="checkbox"/> 20"W x 24"D	ST-SPCGSP2024	£300.00
<input type="checkbox"/> 17"W x 19"D	ST-SPCGSP1719	£250.00	<input type="checkbox"/> 22"W x 18"D	ST-SPCGSP2218	£300.00
<input type="checkbox"/> 18"W x 16"D	ST-SPCGSP1816	£250.00	<input type="checkbox"/> 22"W x 20"D	ST-SPCGSP2220	£300.00
<input type="checkbox"/> 18"W x 17"D	ST-SPCGSP1817	£250.00	<input type="checkbox"/> 22"W x 22"D	ST-SPCGSP2222	£300.00
<input type="checkbox"/> 18"W x 18"D	ST-SPCGSP1818	£250.00	<input type="checkbox"/> 22"W x 24"D	ST-SPCGSP2224	£300.00
<input type="checkbox"/> 18"W x 20"D	ST-SPCGSP1820	£250.00	<input type="checkbox"/> 24"W x 18"D	ST-SPCGSP2418	£300.00
<input type="checkbox"/> 18"W x 22"D	ST-SPCGSP1822	£250.00	<input type="checkbox"/> 24"W x 20"D	ST-SPCGSP2420	£300.00
<input type="checkbox"/> 19"W x 17"D	ST-SPCGSP1917	£250.00	<input type="checkbox"/> 24"W x 22"D	ST-SPCGSP2422	£300.00
<input type="checkbox"/> 19"W x 18"D	ST-SPCGSP1918	£250.00	<input type="checkbox"/> 24"W x 24"D	ST-SPCGSP2424	£300.00
<input type="checkbox"/> 19"W x 19"D	ST-SPCGSP1919	£250.00			

### Add a Solid Base


- |  |               |   |               |
|--|---------------|---|---------------|
| <input type="checkbox"/> ABS Flat Rigidizer      | £35.00 (each) | <input type="checkbox"/> ACM Rigidizer      | £75.00 (each) |
| <i>Not available with TRU-Comfort 2 cushion.</i> |               | <i>For use with manual wheelchair only.</i> |               |

### Additional Fluid Proof Liners

- |   |               |   |               |
|---|---------------|---|---------------|
| <input type="checkbox"/> Fluid Proof Urethane Liner     | £45.00 (each) | <input type="checkbox"/> Fluid Proof Polyurethane Liner | £45.00 (each) |
| <i>Not available with Simplicity and TRU-Comfort 2.</i> |               | <i>Only available with TRU-Comfort 2.</i>               |               |



## Custom TRU-Comfort 2 SPP Cushion

Item	Part	Price
<input type="checkbox"/> Custom Size	ST-TC2SPPCUST	POA
<input type="checkbox"/> Custom Liner	CUS1709085-CUST	POA
<hr/> Detail		
<hr/>		
<hr/>		
<hr/>		
<hr/>		





## Custom Solution SPP Cushion

Item	Part	Price
<input type="checkbox"/> Custom Size	ST-SOLSPPCUST	POA
<input type="checkbox"/> Custom Liner	CUS1709070-CUST	POA
Detail		



## Custom Spectrum SPP Cushion

Item	Part	Price
<input type="checkbox"/> Custom Gel	ST-SPCGPCUST	POA
<input type="checkbox"/> Custom Foam	ST-SPCFSPPCUST	POA
<input type="checkbox"/> Custom Liner	CUS1709070-CUST	POA
*Custom Spectrum Foam and custom liners over 24" will have a higher Price.		
Detail		



### Price

<b>Sub Total:</b>	
<b>Options:</b>	
<b>Total:</b>	

### Special Instructions

---

---

---

---

---

---

---

---

**Quantum Notes:**

- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%. Overall measurements will vary based on seating and accessory selections.