




# Stealth Solution SPP Cushion Order Form Retail

32 Wedgwood Rd, Bicester, OX26 4UL | 01869 324600 | sales@quantumrehab.co.uk

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_  
 Retailer: \_\_\_\_\_  
 Contact: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
 Email: \_\_\_\_\_  
 PO Number: \_\_\_\_\_  
 Marked for: \_\_\_\_\_  
 Ship to Address: \_\_\_\_\_  
 City: \_\_\_\_\_ Postcode: \_\_\_\_\_



This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Prices, specifications, part numbers, and availability are subject to change without notice.

## SOLUTION SPP CUSHION

This cushion offers skin protection and positioning with a Coolcore, 4-way stretch, breathable cover and an urethane liner. Weight capacities are 79.37 kg for 10"-14" width, 158.75 kg for 15"-21" width, and 226.79 kg for 22" width and greater. The XXYY portion of all part numbers on this page directly relate to the dimensions of the cushion, where XX represents the width and YY represents the depth. For example, a 16" wide and 18" deep cushion would be ST-SOLSPP1618. Available sizes are shown on the chart below.

### 1. Cushion Selection

Solution SPP Cushion, 10"-20" Wide .....£180  
 Solution SPP Cushion, 22"-24" Wide ..... £200  
 Cushion Part: **ST-SOLSPPXXYY**  
 Cushion Width: \_\_\_\_\_  
 Cushion Depth: \_\_\_\_\_  
 Complete with Coolcore 4-way stretch breathable cover and Urethane liner.

### 2. Alternative Cushion Cover

The Solution SPP cushion is standard with a Coolcore cover SOLSPP-NBLK\*\*\*\* which can be substituted with a Mesh or Stealth Tek cover (Please select one)

Alternative Mesh Cover .....£FOC  
 Cover Part: **SOLSPP-MBLK XXYY**  
 Cushion Width: \_\_\_\_\_  
 Cushion Depth: \_\_\_\_\_  
 Alternative Stealth Tek Cover .....£FOC  
 Cover Part: **SOLSPP-TBLK XXYY**  
 Cushion Width: \_\_\_\_\_  
 Cushion Depth: \_\_\_\_\_

### 3. Alternative Fluid Proof Liner

The Solution SPP cushion is standard with a Urethane liner URL1709070\*\*\*\* which can be substituted with a Polyurethane Liner

Alternative Polyurethane Liner .....£48  
 Cover Part: **PUL1709030 XXYY**  
 Cushion Width: \_\_\_\_\_  
 Cushion Depth: \_\_\_\_\_

### 4. Add Additional Spare Covers & Liners

Coolcore Cover SOLSPP-NBLK\*\*\*\* .....£70  
 Stealth Tek Cover SOLSPP-TBLK\*\*\*\* .....£50  
 Mesh Cover SOLSPP-MBLK\*\*\*\* .....£44  
 Polyurethane Liner PUL1709030\*\*\*\* .....£48  
 Urethane Liner URL1709070\*\*\*\* .....£30

		Width										
		10"	12"	14"	15"	16"	17"	18"	19"	20"	22"	24"
Depth	10"			-	-	-	-	-	-	-	-	-
	12"				-	-	-	-	-	-	-	-
	14"					-	-	-	-	-	-	-
	15"						-	-	-	-	-	-
	16"							-	-		-	-
	17"								-	-	-	-
	18"									-	-	-
	19"										-	-
	20"											-
	22"											
24"												



Standard viscoelastic foam top with twin cell aqueous viscous gel (heavier gel) insert in ischial well, mid-level contour

Weight Capacities  
 10-14" - 79.37 kg  
 15-21" - 158.75 kg  
 22" and Above - 226.79 kg

Overall Weight  
 2.4 kg

Overall Height  
 3"

10 x 10" up to 24 x 24"  
 Odd sizes in 17" and 19"  
 Custom sizing available

**5. Add Solid Base**

22" width and greater cushions include a flat rigidizer as standard

**1. ABS Flat Rigidizer**

Cushion Size	Part	Cushion Size	Part	Cushion Size	Part
10"W x 10"D.....	DWR1010E188A ..... £34	16"W x 16"D .....	DWR1010E191G ..... £34	19"W x 19"D ....	MEC1804152-1919..... £34
10"W x 12"D.....	DWR1010E188C ..... £34	16"W x 18"D .....	DWR1010E191I..... £34	20"W x 16"D .....	DWR1010E193G ..... £34
12"W x 10"D.....	DWR1010E189A ..... £34	16"W x 20"D.....	DWR1010E191K ..... £34	20"W x 18"D.....	DWR1010E193I..... £34
12"W x 12"D.....	DWR1010E189C ..... £34	17"W x 17"D.....	DWR1010E234H..... £34	20"W x 20"D.....	DWR1010E193K ..... £34
12"W x 14"D.....	DWR1010E189E ..... £34	17"W x 19"D.....	DWR1010E234J ..... £34	20"W x 22"D.....	DWR1010E193M ..... £34
14"W x 12"D.....	DWR1010E190C ..... £34	18"W x 16"D.....	DWR1010E192G ..... £34	22"W x 18"D.....	DWR1010E194I..... £34
14"W x 14"D.....	DWR1010E190E ..... £34	18"W x 18"D.....	DWR1010E192I..... £34	22"W x 20"D.....	DWR1010E194K ..... £34
14"W x 16"D.....	DWR1010E190G ..... £34	18"W x 20"D.....	DWR1010E192K ..... £34	22"W x 22"D.....	DWR1010E194M ..... £34
14"W x 18"D.....	DWR1010E190I..... £34	18"W x 22"D.....	DWR1010E192M ..... £34	24"W x 20"D.....	DWR1010E236K ..... £34
15"W x 18"D.....	DWR1010E232H..... £34	19"W x 17"D.....	MEC1804152-1917..... £34	24"W x 22"D.....	DWR1010E236M ..... £34

**2. ACM Rigidizer**

For use on manual wheelchairs only

Cushion Size	Part	Cushion Size	Part	Cushion Size	Part
15"W x 18"D .....	CUS142467 ..... £56	18"W x 16"D .....	CUS142466 ..... £56	20"W x 18"D.....	CUS142465 ..... £56
16"W x 16"D .....	CUS142467 ..... £56	18"W x 18"D.....	CUS142466 ..... £56	20"W x 19"D.....	CUS142465 ..... £56
16"W x 17"D .....	CUS142467 ..... £56	18"W x 19"D.....	CUS142466 ..... £56	20"W x 20"D.....	CUS142465 ..... £56
16"W x 18"D .....	CUS142467 ..... £56	18"W x 20"D.....	CUS142466 ..... £56		

**SPECIAL ORDER SOLUTION SPP CUSHION**

Please detail your special order cushion request. All special order cushion requests must be approved with the Quantum Custom Department. Please contact Quantum Inside Sales or Custom Department with questions. Extended lead time will apply. Cushion part numbers involved will be the cushion part number listed in previous sections modified to meet special order request stated in details below. Drawings are available below for sketch details.

Special Order Solution SPP Cushion. ....£350

Details

---



---



---




---



---

Top View



Front View



Right Side View

