

Quantum®

32 Wedgwood Road, Bicester,
Oxon, OX26 4UL
Phone: 01869 324600
Fax: 01869 323070
Email: sales@quantumrehab.co.uk
Web: www.quantumrehab.co.uk



Dealer September 2018 (V1.8)



Aspen with TB Flex Seating

Dealer Name: _____ Order Date: _____
 Address: _____

 _____ Post Code: _____
 Phone _____ Email: _____
 PO No: _____ Order Ref: _____
 Delivery Address: _____

 _____ Post Code: _____
 Quantum Specialist: _____

Standard Features

- Seat Width Adjustment 15" - 21"
- Seat Depth Adjustment 15" - 20"
- 2 x 55ah batteries
- Left or Right Joystick Mount
- Max User Weight 21 Stone 6lbs
- 9" x 3" Front Castors, Flat Free
- 14" x 3" Pneumatic Drive Tyres
- Manual Settable Seat Tilt
- 4 Mph (6 Km/h) Motors
- Rear Suspension
- Height Adjustable Removable Armrests
- 70° Swing-Away Legrests
- Angle Adjustable Flip Forward Backrest
- Tension Adjustable Backrest
- 4 Point Tie Down System

Standard = Standard Factory Feature | FOC = Select instead of Standard Factory Feature

Base Model

Aspen with TB Flex Seating - £3,360.00

Drive Wheels

- 14" x 3" Pneumatic, Black Tyres Standard
- 14" x 3" Flat Free, Black Tyres £150.00

Castor Wheels

- 9" x 3" Flat Free, Black Tyres Standard
- 8" x 2" Flat Free, Black Tyres £55.00

Colours



Seat Dimensions

Seat Width

- 15"
- 16"
- 17"
- 18"
- 19"
- 20"
- 21" Only available with tension adjustable backrest

Seat Depth

- 15"
- 16"
- 17"
- 18"
- 19"
- 20"

Please note that the smallest standard seat cushion available is 16" x 16", and the largest is 20" x 20"

Seat Cushions

- | | |
|--|----------|
| <input type="checkbox"/> Sitrite Premier | Standard |
| <input type="checkbox"/> Sitrite Premier Plus | £75.00 |
| <input type="checkbox"/> TRU-Comfort Plus | £80.00 |
| <input type="checkbox"/> 16"W x 16"D <input type="checkbox"/> 16"W x 18"D <input type="checkbox"/> 16"W x 20"D <input type="checkbox"/> 18"W x 16"D <input type="checkbox"/> 18"W x 18"D | |
| <input type="checkbox"/> 18"W x 20"D <input type="checkbox"/> 20"W x 16"D <input type="checkbox"/> 20"W x 18"D <input type="checkbox"/> 20"W x 20"D | |

Please contact your Quantum Specialist for alternative cushions and sizes

Seat Positioning Adjustments

- | | |
|---|----------|
| <input type="checkbox"/> Manual Tilt 0° to 30° with a telescopic tube | Standard |
| <input type="checkbox"/> TB-Flex® Power Tilt 0° to 30°. Full Drive Lock at 20° when combined with Power Recline. Uses 1 power option. | £450.00 |

Seat to Floor Height

- | | | |
|--------------------------------|---|--|
| <input type="checkbox"/> 17.5" | <input type="checkbox"/> 18" Changes Tilt range to -5° to 25° | <input type="checkbox"/> 18.5" Reduces Tilt range to -10° to 20° |
|--------------------------------|---|--|

Backrest Upholstery



- Tension Adjustable Back**
Standard



- Sitrite Premier**
£135.00



- Sitrite Premier Plus**
£170.00



- ADI**
£.....
Part No:
.....

Range of Adjustment: Sitrite Premier and Premier Plus = 24.5" High to 28" High. (Factory Set at 24.5") Tension Adjustable = 19" to 23"

Backrest Positioning Adjustments

- | | |
|--|----------|
| <input type="checkbox"/> Manual Settable / Flip Forward Range of adjustments 95° to 125° | Standard |
| <input type="checkbox"/> Attendant Controlled Manual Recline Range of adjustments 95° to 125°. Not suitable for occupied transportation. | £150.00 |
| <input type="checkbox"/> Powered Recline Range of adjustments 95° to 125°. Uses 1 power option. | £345.00 |

Armrests

- | | |
|--|----------|
| <input type="checkbox"/> TB Flex Removable Drop-in Armrests Range of adjustments 9" to 14" | Standard |
| <input type="checkbox"/> TB Flex Removable Armrests with Quick Release Levers Range of adjustments 9" to 14" | £85.00 |
| <input type="checkbox"/> 2-Post Flip Up Removable Armrests Range of adjustments 9" to 13" | £280.00 |

Armpads

- | | | | |
|---------------------------------------|---------------------------------------|---|---|
| <input type="checkbox"/> 10" Standard | <input type="checkbox"/> 14" Standard | <input type="checkbox"/> 10" Waterfall
<small>Only available with 2-post flip up</small> | <input type="checkbox"/> 14" Waterfall
<small>Only available with 2-post flip up</small> |
|---------------------------------------|---------------------------------------|---|---|

Joystick Options

- | | |
|---|----------|
| <input type="checkbox"/> NE + Controls | Standard |
| <input type="checkbox"/> VR2 Controls | FOC |
| <input type="checkbox"/> Q-logic 2 Expandable Controls 2 power option ready | £350.00 |
| <input type="checkbox"/> Q-logic 2 Expandable Controls 5 power option ready | £550.00 |
| <input type="checkbox"/> R-Net Expandable Controls 1 power option ready | £350.00 |
| <input type="checkbox"/> R-Net Expandable Controls 5 power option ready | £550.00 |

Joystick Mounting Options

- | | | | | |
|---|----------|---|---------|--------------------------------|
| <input type="checkbox"/> Fixed | Standard | <input type="checkbox"/> Retract 4 Swing-away Joystick Arm Kit (Full 180° Movement) | £150.00 | <input type="checkbox"/> Left |
| <input type="checkbox"/> Swing-Away | £75.00 | <input type="checkbox"/> Dual Control | £380.00 | <input type="checkbox"/> Right |
| <input type="checkbox"/> Attendant Control Only | £100.00 | | | |

Legrests

- Seat Pan to Heel Measurement Required**
- | | |
|--|----------|
| <input type="checkbox"/> 70° Swing-away legrests | Standard |
| <input type="checkbox"/> Quantum 60° to 90° Adjustable Swing-away Legrests | £335.00 |
| <input type="checkbox"/> Swing-Away Elevating Legrests | £250.00 |

Headrests

<input type="checkbox"/> Quantum 10" W Headrest	£140.00
<input type="checkbox"/> Stealth Small Contoured 9.75"W x 3"H	£200.00
<input type="checkbox"/> Stealth Large Contoured 11"W x 6"H	£200.00
<input type="checkbox"/> Stealth 10" Comfort Plus	£200.00
<input type="checkbox"/> Stealth 14" Comfort Plus	£200.00
<input type="checkbox"/> Stealth 19" Tri Comfort	£200.00

Please visit www.stealthproducts.com for more headrest options

Posture Supports

Pelvic Belts

<input type="checkbox"/> Auto Buckle Lap Belt, Push Button Release	Standard	
<input type="checkbox"/> iFit 2 Point Padded, Single pull, Medium	<input type="checkbox"/> Push Button Release £60.00 (IFPB5213)	<input type="checkbox"/> Side Release £49.00 (IFPB6213)
<input type="checkbox"/> iFit 2 Point Padded, Single pull, Large	<input type="checkbox"/> Push Button Release £65.00 (IFPB5314)	<input type="checkbox"/> Side Release £55.00 (IFPB6314)
<input type="checkbox"/> iFit 2 Point Padded, Single pull, XLarge	<input type="checkbox"/> Push Button Release £66.00 (IFPB5315)	<input type="checkbox"/> Side Release £55.00 (IFPB6315)

Please refer to iFit catalogue for additional options

Posture Supports (continued)

Stealth Thoracic Lateral Hardware with Backcane Mount and Pad Swing-Away laterals may interfere with flip-up armrest

<input type="checkbox"/> Left: £195.00	<input type="checkbox"/> Curved Pad, 3" x 5"	<input type="checkbox"/> Curved Pad, 4" x 6"
<input type="checkbox"/> Right: £195.00	<input type="checkbox"/> Curved Pad, 3" x 5"	<input type="checkbox"/> Curved Pad, 4" x 6"



Accessories

<input type="checkbox"/> Kerb Climber	£125.00
<input type="checkbox"/> Quick Release Transfer Bars <small>Not Compatible with 2-Post Flip Up Removable armrest</small>	£125.00
<input type="checkbox"/> Quantum Backpack	£40.00
<input type="checkbox"/> Swing-Away Tray Table	£350.00

Price

Sub Total:	£3,360.00
Options:	
Total:	

Notes

1. All measurements given are approximate.
2. By faxing or mailing this order you accept all responsibility for the configuration ordered.
3. Errors made in filling out this order, may result in further charges or delays.
4. If this order has a deadline you must advise us ASAP. We will try to deliver for this date, but cannot accept liability if we are unable to fulfil this order by your deadline.
5. Information relating to crash testing is correct at time of printing.
Pride Mobility constantly reviews and updates this information, if you require further information about crash testing, please contact your Quantum Specialist